

## **U965**Bed Pan & Urine Bottle Rack

	JOB #:
	ITEM #:
	QUANTITY:

Revision 03/11

MOUNTING	HANDING	TYPE	SIZE	OTHER	
	12 1/4 <u>"</u> 311mm	 			1/4"
26" 660mm	4 3/4" 120mm		285mm	5"   127mm   6 3/4"   171mm   20 1/4"   514mm	
		1" (TYP.) 25mm	Mounting Brack	6 3/4" 3 1/4" 171mm 82mm	6" m

## SPECIFICATIONS INSTALLATION

**Back Plate:** Fabricated of #18 ga. stainless steel with a #4 satin finish. **Holding Straps:** Fabricated of #16 ga. stainless steel and finished to match rest of unit.

**Hangers:** Fabricated of #22 ga. stainless steel. **Overall Size:** 12 3/4"W x 26"H x 4 9/16"D

311mm x 660mm x 103mm

 $\textbf{Capacity:} \ \, \textbf{One (1) urine bottle \& one (1) standard bedpan}$ 

Place rack on wall surface making sure unit is level, transfer mounting point locations. Drill pilot holes and secure with proper fasteners per each application. (Fasteners not included)

## NOTES

|